

CURITY^{*} Non-Adhering

Oil Emulsion Dressing



tyco

Healthcare

KENDALL

CURITY* N-A

Product Features

Non-Adherent	Reduces pain and trauma during dressing changes Will not disrupt healing process
Oil Emulsion Impregnate	Non-adhering Allows exudates to readily pass through to absorbent layer without occlusion Helps prevent pooling of fluid and resulting maceration Oil emulsion blend is inert, thus it is compatible with other treatments
Open Mesh, Knitted Fabric Dressing (Consistent, porous construction)	Allow exudate to readily to pass through to absorbent layer without occlusion Will not shred or unravel when trimmed to size
Fewer dressing folds	Prevents uneven distribution of impregnate at the fold which may contribute to maceration of skin

Composition and Characteristic

CURITY Oil Emulsion Dressings are non-adherent, acetate dressings with porous construction, impregnated with oil emulsion blend.

The product is available in strips (tube packing) and dressing configurations.

As a primary wound contact layer, it covers and protects wounds against environmental contamination and allows exudates to readily pass through to an absorbent secondary dressing without occlusion. It may be used in combination with other dressings. Effective coverage is obtained with a single layer of this dressing.

Indications

- Minimal to heavy exuding wounds
- Burns
- Skin grafts/reconstructive surgery
- Abrasions
- Skin tears
- Staple / suture lines

Ordering Information

CURITY* N-A

	REF	Pack Quantity	Case Quantity
1cmx3,6m	6111	1	12
8cmx8cm	6112	50	600
8cmx20cm	6113	24	144
8cmx41cm	6114	36	216
8cmx20cm (3xpkg)	6115	36	216
13cmx23cm	6116	12	72

Test data and study details on file. Reproduced with kind permission of Prof. Elia Ricci, University of Turin, ITALY
* CURITY is a trademark of Tyco Healthcare Group LP. © 2003 Tyco Healthcare Group LP. All rights reserved.

P-WC-CuriNA/GB

tyco
Healthcare

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334

Tyco Healthcare Regional Export
154 Fareham Road
GOSPORT
Hampshire
PO13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare
Birchwood Court
Montrose Road
P.O. Box 8108
Halfway House
1685 MIDRAND
SOUTH AFRICA
Tel: +27 (0)11 265 2000
Fax: +27 (0)11 265 2050

Tyco Healthcare Nederland BV
Hogeweg 105
530 1 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 93
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare Norden AB
Gårdsvägen 18
SE-169 67 SOLNA
SWEDEN
Tel: +46 (0) 8 585 605 05
Fax: +46 (0) 8 585 605 06

www.tycohealthcare.com

KENDALL