

Cardiology & Thoracic Drainage

ARGYLE*

A COMPLETE
RANGE OF
SAFE
SOLUTIONS



tyco

Healthcare

Chest Drainage

- THORA-SEAL I Chest Drainage Unit
- SINGLE SEAL Chest Drainage Unit
- THORA-SEAL II Chest Drainage Unit
- THORA-SEAL III Chest Drainage Unit
- THORA-SEAL III Autotransfusion Chest Drainage Unit
- THORA-SEAL III Replacement Collection Chamber
- THORA-SEAL III Autotransfusion Chest Drainage System
- ARGYLE AQUA-SEAL Chest Drainage Unit
- ARGYLE AQUA-SEAL Dual Drain Chest Drainage Unit
- ARGYLE AQUA-SEAL Autotransfusion Chest Drainage Unit & System
- ARGYLE AQUA-SEAL All Purpose Autotransfusion Chest Drainage Unit
- ARGYLE Autotransfusion Accessory Unit
- ARGYLE AQUA-SEAL Evacuating Autotransfusion Accessory
- ARGYLE SENTINEL SEAL Chest Drainage Unit
- ARGYLE SENTINEL SEAL Dual Drain Chest Drainage Unit
- ARGYLE SENTINEL SEAL Autotransfusion Chest Drainage Unit & System

- ARGYLE DOUBLE SEAL Chest Drainage Unit
- CLINICAL PRODUCTS Bottle Drainage Systems & Accessories

Thoracic Catheters

- ARGYLE Thoracic Catheters
- ARGYLE Silicone Thoracic Catheters
- ARGYLE Trocar Catheters

- ARGYLE Silicone Mediastinal Catheters
- ARGYLE Tubing Extension Sets

- ARGYLE Trocar Catheter Kits and Thoracostomy Procedure Tray

- ARGYLE TURKEL Thoracentesis Set

- ARGYLE Left Ventricular Sump Vent Catheters
- ARGYLE Ferguson Left Vent Catheters

Thoracic Catheter Accessories

- ARGYLE THI Aortic Perfusion Cannulas
- ARGYLE THI Perfusion Adapters
- ARGYLE Parallel “Y” Connector
- MINIDRAIN Pumps
- ARGYLE Gott Aortic Aneurysm Shunts
- ARGYLE Carotid Artery Shunt Kits
- ARGYLE Hepatic Venous Shunt
- ARGYLE Vascular Tourniquet Kit
- ARGYLE SENTINEL LOOPS Silicone Vascular Tapes
- ARGYLE CARDIOPUNCH Disposable Aortic Punch
- ARGYLE Penrose Tubing
- ARGYLE SARATOGA Sump Drains
- ARGYLE MEDICUT Intravenous Cannulas
- ARGYLE EZ-FLO Stopcocks

AQUA-SEAL, ARGYLE, CLINICAL PRODUCTS, DOUBLE SEAL, EZ-FLO, MEDICUT, MINIDRAIN, MONOJECT, SARATOGA, SENTINEL EYE, SENTINEL LINE, SENTINEL LOOPS, SENTINEL SEAL, SINGLE SEAL, THORA-SEAL and TURKEL are registered trademarks of Sherwood Services AG. CARDIOPUNCH is a registered trademark of American Cyanamid Company. TRU-CLOSE and VASCU-FLO are registered trademarks of UreSil, L.P. VELCRO is a registered trademark of Velcro Industries B.V.



THORA-SEAL I Chest Drainage Unit

- Disposable, 2000ml single bottle unit.
- Can be used with gravity or regulated suction.
- 150cm clear PVC tubing connects to thoracic catheter.
- Self-contained floor stand.
- Sterile.

| Code 8884- | Description | Ship Case |
|------------|----------------------------------|-----------|
| 713100 | THORA-SEAL I Chest Drainage Unit | 6 |

| Code 1814- | Description | Ship Case |
|------------|-------------------------------------|-----------|
| 713105 | THORA-SEAL I Bottle only (With Cap) | 6 |
| 713127 | Carrying Frame for 2 bottles | 1 |



SINGLE SEAL Chest Drainage Unit

- Compact, two bottle chest drainage system.
- With Kraton tube and connector.
- Can be used with gravity or regulated suction.
- Specially treated collar reduces foaming.
- Sterile.

| Code 8888- | Description | Graduations | Ship Case |
|------------|------------------------|-------------|-----------|
| 571588 | 2000ml Two Bottle Unit | 50ml | 10 |



THORA-SEAL II Chest Drainage Unit

- Disposable, 2-chamber chest drainage unit.
- Compact, clear break-resistant plastic.
- For a wide range of clinical applications including those with heavy blood loss.
- Non-latex patient tube.
- Self-monitoring 120ml underwater seal chamber prevents progressive resistance.
- Easy-to-read graphics.
- 2600ml capacity.
- Pre-attached floor stand and hangers.
- Sterile.

| Code 1180- | Description | Ship Case |
|------------|-----------------------------------|-----------|
| 571570 | THORA-SEAL II Chest Drainage Unit | 5 |

Chest
Drainage

Chest Drainage



THORA-SEAL III Chest Drainage Unit

- Compact, clear 3-chamber unit.
- Removable collection chamber.
- Non-latex patient tube.
- Automatic positive pressure and negative pressure relief valve.
- Large fill openings for quick set-up.
- Easy-to-read graphics.
- Set-up instructions printed on unit.
- Unique baffle system to prevent mixing of fluids.
- Pre-attached floor stand and hangers.
- Sterile.

| Code 8884- | Description | Ship Case |
|---------------|------------------------------------|--------------|
| 713308 | THORA-SEAL III Chest Drainage Unit | 4 |



THORA-SEAL III Autotransfusion Chest Drainage Unit

- Compact, clear 3-chamber autotransfusion-ready chest drainage unit.
- Compatible with ARGYLE Autotransfusion Accessory Unit.
- Removable collection chamber.
- Non-latex patient tube.
- Automatic positive pressure and negative pressure relief valve.
- Large fill openings for quick set-up.
- Unique baffle system to prevent mixing of fluids.
- Pre-attached floor stand and hangers.
- Colour-coded locking connectors.
- Simple hook and VELCRO design.
- Sterile.

| Code 8884- | Description | Ship Case |
|---------------|--|--------------|
| 713176 | THORA-SEAL III Autotransfusion Chest Drainage Unit | 4 |



THORA-SEAL III Replacement Collection Chamber

- For THORA-SEAL III Chest Drainage and Autotransfusion Units.
- Twist-on cap included.
- Sterile.

| Code 8884- | Description | Ship Case |
|---------------|--------------------------------|--------------|
| 713900 | Replacement Collection Chamber | 6 |



THORA-SEAL III Autotransfusion Chest Drainage System

- Complete pre-attached autotransfusion system.
- Compatible with ARGYLE Autotransfusion Accessory Unit.
- Compact, clear 3-chamber autotransfusion chest drainage unit.
- Removable collection chamber.
- Automatic positive pressure and negative pressure relief valve.
- Large fill openings for quick set-up.
- Easy-to-read graphics.
- Set-up instructions printed on unit.
- Unique baffle system to prevent mixing of fluids.
- Pre-attached floor stand and hangers.
- Non-latex patient tube.
- Colour-coded locking connectors.
- Simple hook and VELCRO design.
- Sterile.

| Code 8884- | Description | Ship Case |
|---------------|---|--------------|
| 713192 | THORA-SEAL III Autotransfusion Chest Drainage Unit with Pre-attached Autotransfusion Accessory Unit | 3 |

Chest
Drainage



ARGYLE AQUA-SEAL Chest Drainage Unit

- Compact, clear, 3-chamber chest drainage unit.
- Non-latex patient tube.
- Fluid level indicators identify drainage volume.
- Wide, rotational floor stand and hangers.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured water seal.
- Set-up instructions printed on unit.
- Sterile.



ARGYLE AQUA-SEAL Dual Drain Chest Drainage Unit

- Two independent collection chambers and two collection tubes.
- Non-latex patient tubes.
- Compact, clear, 3-chamber chest drainage unit.
- Fluid level indicators identify drainage volume.
- Wide, rotational floor stand and hangers.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured water seal.
- Set-up instructions printed on unit.
- Sterile.

Chest Drainage

| Code 8888- | Description | Ship Case |
|------------|-------------------------------|-----------|
| 571299 | AQUA-SEAL Chest Drainage Unit | 5 |

| Code 8888- | Description | Ship Case |
|------------|--|-----------|
| 571406 | AQUA-SEAL Dual Drain Chest Drainage Unit | 5 |



ARGYLE AQUA-SEAL Autotransfusion Chest Drainage Unit

- Compact, clear, 3-chamber autotransfusion-ready chest drainage unit.
- Compatible with ARGYLE Autotransfusion Accessory Unit.
- Non-latex patient tube.
- Fluid level indicators identify drainage volume.
- Wide, rotational floor stand and hangers.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured water seal.
- Set-up instructions printed on unit.
- Colour-coded locking connectors.
- Simple hook and VELCRO design.
- Sterile.

| Code 8888- | Description | Ship Case |
|---------------|---|--------------|
| 571315 | AQUA-SEAL Autotransfusion Chest Drainage Unit | 5 |



ARGYLE AQUA-SEAL Autotransfusion Chest Drainage System

- Complete pre-attached autotransfusion system.
- Compatible with ARGYLE Autotransfusion Accessory Unit.
- Compact, clear, 3-chamber autotransfusion chest drainage unit.
- Non-latex patient tube.
- Fluid level indicators identify drainage volume.
- Wide, rotational floor stand and hangers.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured water seal.
- Set-up instructions printed on unit.
- Colour-coded locking connectors.
- Simple hook and VELCRO design.
- Sterile.

| Code 8888- | Description | Ship Case |
|---------------|---|--------------|
| 571323 | AQUA-SEAL Autotransfusion Chest Drainage Unit w/Pre-attached Autotransfusion Accessory Unit | 3 |

Chest
Drainage

Chest Drainage



ARGYLE AQUA-SEAL All Purpose Autotransfusion Chest Drainage Unit

- Three convenient post-operative ATS options including continuous, post-collection batch evacuation and in-line collection.
- Out-the-top ATS Infusion Line for clean, convenient blood retrieval.
- Exclusive Redi-Prime Valve for access to ATS Infusion Line.
- Clot Chamber for clot volume assessment.
- Non-latex patient tube and components.
- Automatic/Manual Negative Pressure Relief.
- Blue-coloured water seal.
- Wide rotational floorstand and hangers.
- Flow control valve for convenient suction regulation.
- Sterile.

| Code 8888- | Description | Ship Case |
|------------|--|-----------|
| 571364 | AQUA-SEAL All Purpose Auto-transfusion Chest Drainage Unit | 5 |



ARGYLE Autotransfusion Accessory Unit

- Compatible with THORA-SEAL III Autotransfusion Chest Drainage Unit, THORA-SEAL III Autotransfusion Chest Drainage System, AQUA-SEAL Autotransfusion Chest Drainage Unit and AQUA-SEAL Autotransfusion Chest Drainage System.
- Fits easily into standard pressure infusion cuff.
- Automatic blood overflow to chest drainage unit.
- Graduations printed on the autotransfusion bag.
- Can be used with gravity or continuous reinfusion.
- Injection port for addition of anti-coagulants.
- Non-latex tubing.
- Sterile.

| Code 8884- | Description | Ship Case |
|------------|--|-----------|
| 713184 | ARGYLE Autotransfusion receptacle for collection and reinfusion of blood | 5 |

ARGYLE AQUA-SEAL Evacuating Autotransfusion Accessory

- Exclusively used with AQUA-SEAL All Purpose Autotransfusion Chest Drainage Unit.
- Quick, simple post collection batch evacuation for autotransfusion.
- Vacuum-controlled filling process.
- Blood bag feels and handles like conventional blood bag.
- Optional Refill Bag.
- Sterile.

| Code 8884- | Description | Ship Case |
|------------|--|-----------|
| 571368 | AQUA-SEAL Evacuating Auto-transfusion Refill Bag | 5 |



ARGYLE SENTINEL SEAL Chest Drainage Unit

- Compact, clear, 4-chamber chest drainage unit.
- Latex-free tubing.
- Precision suction regulator.
- Exclusive patient assessment chamber monitors actual negative intrathoracic pressure.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured patient assessment chamber.
- Set-up instructions printed on unit.
- Retractable floor stand and hangers.
- Sterile.

| Code 8888- | Description | Ship Case |
|---------------|-----------------------------------|--------------|
| 571562 | SENTINEL SEAL Chest Drainage Unit | 5 |



ARGYLE SENTINEL SEAL Dual Drain Chest Drainage Unit

- Two independent collection chambers and two Kraton collection tubes.
- Compact, clear, 4-chamber chest drainage unit.
- Precision suction regulator.
- Exclusive patient assessment chamber monitors actual negative intrathoracic pressure.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured patient assessment chamber.
- Set-up instructions printed on unit.
- Retractable floor stand and hangers.
- Sterile.

| Code 8888- | Description | Ship Case |
|---------------|--|--------------|
| 571513 | SENTINEL SEAL Dual Drain Chest Drainage Unit | 5 |

Chest
Drainage

Chest Drainage



ARGYLE SENTINEL SEAL Autotransfusion Chest Drainage Unit

- Compact, clear, 4-chamber autotransfusion-ready chest drainage unit.
- Precision suction regulator.
- Exclusive patient assessment chamber monitors actual negative intrathoracic pressure.
- Clear, easy-to-read graphics.
- Paediatric and adult graduations.
- Automatic positive pressure relief valve.
- Filtered negative pressure relief vent.
- Blue-coloured patient assessment chamber.
- Set-up instructions printed on unit.
- Retractable floor stand and hangers.
- Colour-coded locking connectors.
- Quick release latch design.
- Compatible with SENTINEL SEAL Autotransfusion Accessory Unit.
- Sterile.

| Code 8888- | Description | Ship Case |
|------------|---|-----------|
| 571489 | SENTINEL SEAL Autotransfusion Chest Drainage Unit | 5 |



ARGYLE SENTINEL SEAL Autotransfusion Chest Drainage System

- Complete pre-attached autotransfusion system.
- Compact, clear, 4-chamber autotransfusion chest drainage unit.
- Compatible with SENTINEL SEAL Autotransfusion Accessory Unit.
- Precision suction regulator.
- Exclusive patient assessment chamber monitors actual negative intrathoracic pressure.
- Easy-to-read graphics.
- Paediatric and adult graduations.
- Filtered negative pressure relief vent.
- Blue-coloured water.
- Set-up instructions printed on unit.
- Retractable floor stand and hangers.
- Colour-coded locking connectors.
- Quick-release latch design.
- Sterile.

| Code 8888- | Description | Ship Case |
|------------|---|-----------|
| 571463 | SENTINEL SEAL Autotransfusion Chest Drainage Unit w/Pre-attached Autotransfusion Accessory Unit | 3 |



**ARGYLE SENTINEL SEAL
Autotransfusion Accessory Unit**

— Use Code 8884713184



**ARGYLE DOUBLE SEAL Chest
Drainage Unit**

- Exclusive patient assessment chamber monitors actual negative intrathoracic pressure.
- Automatic negative pressure relief system.
- Blue-coloured water.
- Includes floor stand and hangers.
- Sterile.

| Code 8888- | Description | Ship Case |
|---------------|------------------------------------|--------------|
| 571547 | DOUBLE SEAL Chest Drainage Unit | 5 |

Chest
Drainage

Chest Drainage



CLINICAL PRODUCTS One and Two Bottle Chest Drainage Systems

- Crystal clear glass bottles.
- Accurate graduations.
- 150cm clear PVC tubing connects to thoracic catheter.
- Sterile.

| Code 8884- | Description | Graduations | Ship Case |
|---------------------------|--------------------------------------|-------------|-----------|
| 722100 | 2000ml Bottle only | 50ml | 6 |
| One Bottle Systems | | | |
| 726100 | Standard 2000ml Bottle | 50ml | 6 |
| 726500 | Dual Drainage 2000ml Bottle | 50ml | 6 |
| 728100 | Infant 90ml Bottle with VELCRO strap | 1ml | 12 |
| Two Bottle Systems | | | |
| 726200 | Standard 2000ml Bottles | 50ml | 4 |
| 726400 | Dual Drainage 2000ml Bottles | 50ml | 4 |



CLINICAL PRODUCTS Glass Bottle Chest Drainage Tubing Sets

- For any size glass bottle with 38mm neck.
- Sterile.

| Code 8884- | Description | Ship Case |
|------------|----------------------------------|-----------|
| 729300 | Fixed One Bottle Tubing Set | 12 |
| 729500 | Adjustable One Bottle Tubing Set | 12 |
| 729600 | Adjustable Two Bottle Tubing Set | 12 |



CLINICAL PRODUCTS Floor Rack for Glass Bottles

- Solid wire construction.
- Lightweight.

| Code 8884- | Description | Ship Case |
|------------|-------------|-----------|
| 722400 | Floor Rack | 12 |



ARGYLE Thoracic Catheters

- Clear thermosensitive PVC catheter.
- Radiopaque SENTINEL LINE with SENTINEL EYE.
- Smooth side eyes and open distal tip.
- Numeric depth marks at 2cm increments.
- Bevelled connector end can be cropped to form integral funnel connector.
- Sterile.

| Code 8888- | French Size | O.D. mm | Side Eyes | Length | Ship Case |
|--------------------|----------------|------------|--------------|--------|--------------|
| Straight | | | | | |
| 570507 | 12 | 4.0 | 4 | 51cm | 10 |
| 570515 | 16 | 5.3 | 4 | 51cm | 10 |
| 570523 | 20 | 6.7 | 5 | 51cm | 10 |
| 570531 | 24 | 8.0 | 6 | 51cm | 10 |
| 570549 | 28 | 9.3 | 6 | 51cm | 10 |
| 570556 | 32 | 10.7 | 6 | 51cm | 10 |
| 570564 | 36 | 12.0 | 6 | 51cm | 10 |
| 570572 | 40 | 13.3 | 6 | 51cm | 10 |
| Right Angle | | | | | |
| 571018 | 16 | 5.3 | 4 | 51cm | 10 |
| 571026 | 20 | 6.7 | 5 | 51cm | 10 |
| 571034 | 24 | 8.0 | 6 | 51cm | 10 |
| 571042 | 28 | 9.3 | 6 | 51cm | 10 |
| 571059 | 32 | 10.7 | 6 | 51cm | 10 |
| 571067 | 36 | 12.0 | 6 | 51cm | 10 |
| 571075 | 40 | 13.3 | 6 | 51cm | 10 |



ARGYLE Silicone Thoracic Catheters

- Biocompatible, thrombo-resistant silicone catheter.
- Soft, kink-resistant.
- Radiopaque SENTINEL LINE with SENTINEL EYE.
- Smooth side eyes and open distal tip.
- Bevelled end.
- Sterile.

| Code 8888- | French Size | O.D. mm | Side Eyes | Length | Ship Case |
|--------------------|----------------|------------|--------------|--------|--------------|
| Straight | | | | | |
| 572537 | 20 | 6.7 | 5 | 51cm | 10 |
| 572545 | 24 | 8.0 | 6 | 51cm | 10 |
| 572552 | 28 | 9.3 | 6 | 51cm | 10 |
| 572560 | 32 | 10.7 | 6 | 51cm | 10 |
| 572578 | 36 | 12.0 | 6 | 51cm | 10 |
| 572586 | 40 | 13.3 | 6 | 51cm | 10 |
| Right Angle | | | | | |
| 573030 | 24 | 8.0 | 6 | 51cm | 10 |
| 573048 | 28 | 9.3 | 6 | 51cm | 10 |
| 573055 | 32 | 10.7 | 6 | 51cm | 10 |
| 573063 | 36 | 12.0 | 6 | 51cm | 10 |

Thoracic
Catheters



ARGYLE Trocar Catheters

- Clear thermosensitive PVC catheter over precision aluminum trocar.
- Surgically sharp trocar tip.
- Radiopaque SENTINEL LINE with SENTINEL EYE.
- Numeric depth markings at 2cm increments.
- Integral funnel end.
- Double lumen has one-way injection port.
- Sterile.

| Code 8888- | French Size | O.D. mm | Length | Ship Case |
|------------|-------------|---------|--------|-----------|
| 560805 | 8 | 2.7 | 23cm | 10 |
| 561019 | 10 | 3.3 | 23cm | 10 |
| 561027 | 12 | 4.0 | 23cm | 10 |
| 561035 | 16 | 5.3 | 25cm | 10 |
| 561043 | 20 | 6.7 | 41cm | 10 |
| 561050 | 24 | 8.0 | 41cm | 10 |
| 561068 | 28 | 9.3 | 41cm | 10 |
| 561076 | 32 | 10.7 | 41cm | 10 |

Trocar Replacement Catheters

| | | | | |
|--------|----|------|------|----|
| 561514 | 10 | 3.3 | 23cm | 10 |
| 561522 | 12 | 4.0 | 23cm | 10 |
| 561530 | 16 | 5.3 | 23cm | 10 |
| 561548 | 20 | 6.7 | 25cm | 10 |
| 561555 | 24 | 8.0 | 41cm | 10 |
| 561563 | 28 | 9.3 | 41cm | 10 |
| 561571 | 32 | 10.7 | 41cm | 10 |

| Code 1180- | French Size | O.D. mm | Length | Ship Case |
|------------|-------------|---------|--------|-----------|
|------------|-------------|---------|--------|-----------|

Double Lumen Thoracic Trocar Catheters


| | | | | |
|--------|----|-----|------|----|
| 562009 | 12 | 4.0 | 23cm | 10 |
| 562017 | 16 | 5.3 | 23cm | 10 |
| 562025 | 20 | 6.7 | 25cm | 10 |
| 562033 | 24 | 8.0 | 41cm | 10 |
| 562041 | 28 | 9.3 | 41cm | 10 |

| 1180- | Size | mm | Length | Case |
|-------|------|----|--------|------|
|-------|------|----|--------|------|

Round Tip Trocar Catheters

| | | | | |
|--------|----|-----|------|----|
| 561415 | 10 | 3.3 | 23cm | 10 |
| 561423 | 12 | 4.0 | 23cm | 10 |
| 561449 | 16 | 5.3 | 25cm | 10 |
| 561464 | 20 | 6.7 | 41cm | 10 |

Other Sizes Available On Request



ARGYLE Silicone Mediastinal Catheters

- Biocompatible, thrombo-resistant silicone catheter.
- Flat drainage section with six eyes on one side.
- Round section exits wound.
- Bevelled end.
- Sterile.

| Code 8888- | French Size | O.D. mm | Length | Flat Drain Width | Ship Case |
|------------|-------------|---------|--------|------------------|-----------|
| 573550 | 28 | 9.3 | 51cm | 11.0mm | 10 |
| 573568 | 34 | 11.3 | 51cm | 13.4mm | 10 |

Thoracic Catheters



ARGYLE Tubing Extension Sets

- 90cm Kraton tubing.
- Croppable step connectors.
- Sterile.

| Code 8884- | Description | Ship Case |
|------------|---|-----------|
| 714100 | For use with a single thoracic catheter | 24 |
| 714200 | For use with dual thoracic catheters | 24 |



ARGYLE Trocar Catheter Kits and Thoracostomy Procedure Tray

Sterile Contents:

Prep Side

- Pair Non-Latex Gloves.
- Prep Swabs x 2.
- Prep Cups x 2.
- Fenestrated Drape with Self-Adhesive Attachment Tabs (55cm x 66cm).
- Roll Tape (5cm x 90cm).
- Medicine Cup (50ml).

Procedure Side

- MONOJECT 12ml Syringe with Caps x 2.
- MONOJECT 25ga x 5/8" Needle.
- MONOJECT 20ga x 3 1/2" Needle.
- MONOJECT 20ga x 1 1/2" Needle x 2.
- Scalpel with #11 Blade.
- 2 1/2" #96 Keith Cutting Needle with Attached 0 Silk Suture.
- Gauze Sponges (10cm x 10cm) x 4.
- Pre-cut Gauze Dressings (10cm x 10cm) x 2.
- ARGYLE Trocar Catheter.

| Code 8888- | Catheter French Size | Catheter O.D. mm | Catheter Length | Ship Case |
|---------------|-------------------------|---------------------|--------------------|--------------|
| 565028 | 12 | 4.0 | 23cm | 10 |
| 565036 | 20 | 6.7 | 25cm | 10 |

Thoracic
Catheters

Thoracic Catheters



ARGYLE TURKEL Thoracentesis Set

- Minimises risk of potential pneumothorax.
- One-way valve in needle hub allows air and fluid to escape from pleural space but prevents atmospheric air from entering.
- Safety cannula housed within sharp, bevelled 16ga hollow needle – withdraws into needle shaft during penetration of chest wall.
- Colour-coded change indicator provides visual confirmation of cannula position.
- Sterile contents:
 - Thoracentesis Needle.
 - Polyurethane Cannula.
 - 2000ml Drain Bag.
 - 50ml Luer Lock Syringe.

| Code 8888- | Catheter French Size | Catheter O.D. mm | Ship Case |
|---------------|-------------------------|---------------------|--------------|
| 566067 | 8 | 2.7 | 10 |

Temporarily Discontinued Re-Launch 01.03.05



ARGYLE Left Ventricular Sump Vent Catheters

- Double lumen sump tube.
- Clear PVC catheter.
- Bullet tip.
- With Cath-Lok retention rings.
- X-ray opaque SENTINEL LINE.
- Four side eyes.
- Integral tapered connector end.
- Blue pigtail funnel connector on secondary lumen.
- Sterile.

| Code 8888- | French Size | O.D. mm | Length | Ship Case |
|---------------|----------------|------------|--------|--------------|
| 590208 | 14 | 4.7 | 51cm | 10 |
| 590216 | 18 | 6.0 | 51cm | 10 |
| 590224 | 20 | 6.7 | 51cm | 10 |



ARGYLE Ferguson Left Vent Catheters

- Malleable wire insert.
- Clear PVC catheter.
- Bullet tip.
- X-ray opaque SENTINEL LINE.
- Four side eyes.
- Integral tapered connector end.
- Sterile.

| Code 8888- | French Size | O.D. mm | Length | Ship Case |
|---------------|----------------|------------|--------|--------------|
| 590018 | 18 | 6.0 | 51cm | 10 |
| 590026 | 24 | 8.0 | 51cm | 10 |

Thoracic
Catheters



ARGYLE THI Aortic Perfusion Cannulas

- Clear PVC catheter.
- Blue orientation line runs the length of the cannula.
- Smooth polished bevelled tip.
- Integral funnel connector.
- 51cm cannula length.
- With Cath-Lok retention ring.
- Sterile.

| Code 8888- | Fr. Size | O.D. mm | Length of Tip Beyond Ring | Proximal Lumen | Ship Case* |
|------------|----------|---------|---------------------------|----------------|------------|
|------------|----------|---------|---------------------------|----------------|------------|

Angled Tip

| | | | | | |
|--------|----|-----|------|-------|-----|
| 591008 | 10 | 3.3 | 19mm | 6.4mm | 100 |
| 591123 | 12 | 4.0 | 20mm | 6.4mm | 100 |
| 591024 | 14 | 4.7 | 22mm | 6.4mm | 100 |
| 591149 | 16 | 5.3 | 23mm | 6.4mm | 100 |
| 591040 | 18 | 6.0 | 25mm | 6.4mm | 100 |
| 591065 | 21 | 7.0 | 31mm | 9.5mm | 100 |
| 591081 | 24 | 8.0 | 38mm | 9.5mm | 100 |

Straight Tip

| | | | | | |
|--------|----|-----|------|-------|-----|
| 591016 | 10 | 3.3 | 13mm | 6.4mm | 100 |
| 591131 | 12 | 4.0 | 14mm | 6.4mm | 100 |
| 591032 | 14 | 4.7 | 15mm | 6.4mm | 100 |
| 591156 | 16 | 5.3 | 16mm | 6.4mm | 100 |
| 591057 | 18 | 6.0 | 17mm | 6.4mm | 100 |
| 591073 | 21 | 7.0 | 22mm | 9.5mm | 100 |
| 591099 | 24 | 8.0 | 25mm | 9.5mm | 100 |

*10 Boxes of 10.



ARGYLE THI Perfusion Adapters

- Polypropylene construction.
- Standard barbed male connector with male luer for 9.5mm or 6.4mm connection.
- Sterile.

| Code 8888- | Type | Fits Tubing I.D. | Ship Case* |
|------------|------------|------------------|------------|
| 593012 | Adult | 9.5mm | 100 |
| 593038 | Paediatric | 6.4mm | 100 |

*10 Boxes of 10



ARGYLE Parallel "Y" Connector

- For use with Vena Caval Catheters.
- Durable, clear polycarbonate construction.
- 22mm parallel leg separation.
- Smooth, contoured interior walls.
- Ribbed connector ends.
- Sterile.

| Code 8888- | Size | Ship Case |
|------------|------------------|-----------|
| 584912 | 9.5 x 9.5 x 13mm | 5 |

Thoracic Catheter Accessories

MINIDRAIN Pumps

| Code | Description | Ship Case |
|--------|------------------|-----------|
| 571548 | MINIDRAIN I (6L) | 1 |

Thoracic Catheter Accessories



ARGYLE Gott Aortic Aneurysm Shunts

- Precision-extruded medical grade, clear PVC shunt.
- TDMAC/Heparin-coated* inner lumen and outer surface of catheter tips.
- Open, bevel cut proximal end with smooth side eye.
- Distal end tapers to bevel cut cannula.
- Smooth tapered inner surface.
- Sterile.

| Code 8888- | Tip O.D. | Length | Ship Case |
|------------|----------|--------|-----------|
| 551010 | 7mm | 66cm | 1 |
| 551036 | 9mm | 90cm | 1 |

* *Tridodecylmethylammonium Chloride chemically bonded to Heparin. A process developed at Battelle Memorial Laboratories.*



ARGYLE Hepatic Venous Shunt

- For use in liver transplants.
- Consists of three 9.5mm I.D. catheters joined by a special "Y" connector.
- TDMAC/Heparin-coated* inner lumen of shunt and outer 20cm of the three cannula tips.
- Precision-extruded medical grade, crystal clear PVC shunt.
- Sterile.

| Code 8888- | Size | Ship Case |
|------------|-------|-----------|
| 552018 | Adult | 1 |

* *Tridodecylmethylammonium Chloride chemically bonded to Heparin. A process developed at Battelle Memorial Laboratories.*



ARGYLE Carotid Artery Shunt Kits

- Four French sizes in each package.
- Precision-extruded medical grade, clear PVC shunt.
- X-ray opaque SENTINEL LINE.
- Polished tapered ends.
- Sterile.

| Code 8888- | Contents | Length | Ship Case |
|------------|---|--------|-----------|
| 577775 | 10 Fr (3.3mm) 12 Fr (4.0mm) 14 Fr (4.7mm) | 15cm | 10 |



ARGYLE Vascular Tourniquet Kit

- Includes Red Rubber Sheath with Retaining Plug, an Obturator and Umbilical Tape.
- Used to secure Umbilical Tape, Suture, or SENTINEL LOOPS Tapes.
- Sterile.

| Code 8888- | Contents | Size | Ship Case |
|------------|---|------|-----------|
| 585000 | 2 Vascular Tourniquet Sheaths with plugs and obturators and 2 cut lengths of umbilical tape (3mm x 260cm) | 13cm | 10 |



ARGYLE SENTINEL LOOPS Silicone Vascular Tapes

- Soft, silicone formulation.
- Radiopaque.
- Colour-coded.
- Elliptical cross section.
- Smooth surface.
- Sterile.

| Code 8888- | Colour | Cross Section | Length | Ship Case* |
|---------------------|--------|---------------|--------|------------|
| Regular Size | | | | |
| 585208 | Blue | 1.0 x 2.5mm | 38.0cm | 100 |
| 585216 | Red | 1.0 x 2.5mm | 38.0cm | 100 |
| 585224 | Yellow | 1.0 x 2.5mm | 38.0cm | 100 |
| 585232 | White | 1.0 x 2.5mm | 38.0cm | 100 |
| Small Size | | | | |
| 585240 | Blue | 0.8 x 1.5mm | 35.5cm | 100 |
| 585257 | Red | 0.8 x 1.5mm | 35.5cm | 100 |
| 585265 | Yellow | 0.8 x 1.5mm | 35.5cm | 100 |
| 585273 | White | 0.8 x 1.5mm | 35.5cm | 100 |

*10 Boxes of 10


Thoracic
Catheter
Accessories



ARGYLE CARDIOPUNCH Disposable Aortic Punch

- Ergonomic, hand-contoured pistol grip.
- Cuts even calcified aortic tissue.
- Long barrel provides clear field of vision.
- Produces numerous clean cuts in patients requiring more than one bypass.
- Sterile.

| Code 8886- | Punch Size Diameter | Ship Case |
|---------------|------------------------|--------------|
| 822726 | 2.6mm | 6 |
| 822740 | 4.0mm | 6 |
| 822745 | 4.5mm | 6 |
| 822748 | 4.8mm | 6 |
| 822750 | 5.0mm | 6 |



ARGYLE Penrose Tubing

- Latex construction.
- Radiopaque.
- Sterile.

| Code 8888- | Flat Size | Diameter | Length | Ship Case |
|---------------|--------------|----------|--------|--------------|
| 513002 | 11mm | 6mm | 30cm | 50 |
| 513200 | 13mm | 8mm | 30cm | 50 |
| 513408 | 16mm | 10mm | 30cm | 50 |
| 513606 | 19mm | 13mm | 30cm | 50 |
| 513804 | 25mm | 16mm | 30cm | 50 |
| 514000 | 32mm | 19mm | 30cm | 50 |
| 514208 | 35mm | 22mm | 30cm | 50 |
| 514406 | 38mm | 25mm | 30cm | 50 |
| 514604 | 11mm | 6mm | 46cm | 50 |
| 514802 | 13mm | 8mm | 46cm | 50 |
| 515007 | 16mm | 10mm | 46cm | 50 |
| 515205 | 19mm | 13mm | 46cm | 50 |
| 515403 | 25mm | 16mm | 46cm | 50 |
| 515601 | 32mm | 19mm | 46cm | 50 |
| 515809 | 38mm | 25mm | 46cm | 50 |

Thoracic Catheter Accessories



ARGYLE SARATOGA Sump Drains

- Clear thermosensitive PVC catheter.
- Two sump vents.
- Radiopaque SENTINEL LINE with SENTINEL EYE.
- Bevelled connector end can be cropped to form integral funnel connector.
- Sterile.

| Code 8888- | French Size | O.D. mm | Eye Distribution | Length | Ship Case |
|------------|-------------|---------|------------------|--------|-----------|
| 530519 | 16 | 5.3 | | 30cm | 10 |
| 530527 | 20 | 6.7 | From | 30cm | 10 |
| 530543 | 24 | 8.0 | Distal Tip | 30cm | 10 |
| 530535 | 28 | 9.3 | to 5cm | 30cm | 10 |
| 530550 | 34 | 11.3 | | 30cm | 10 |
| 530568 | 40 | 13.3 | | 30cm | 10 |
| 531012 | 16 | 5.3 | | 51cm | 10 |
| 531020 | 20 | 6.7 | From | 51cm | 10 |
| 531038 | 24 | 8.0 | Distal Tip | 51cm | 10 |
| 531046 | 28 | 9.3 | to 10cm | 51cm | 10 |
| 531053 | 34 | 11.3 | | 51cm | 10 |
| 531061 | 40 | 13.3 | | 51cm | 10 |



ARGYLE MEDICUT Intravenous Cannulas

- One-piece cannula with integral tapered end.
- Thin wall polypropylene construction with smooth interior taper.
- Drawn distal tip.
- 21ga/2ml MONOJECT syringe and stainless steel needle.
- Sterile.

| Code 8888- | Cannula Gauge | Cannula Length | Colour Code | Ship Case |
|------------|---------------|----------------|-------------|-----------|
| 100206 | 12 (2.8mm) | 5.0cm | Light Blue | 50 |
| 100214 | 14 (2.1mm) | 5.0cm | Tan | 50 |
| 100222 | 16 (1.7mm) | 5.0cm | Grey | 50 |
| 100230 | 18 (1.2mm) | 5.0cm | Green | 50 |
| 100248 | 20 (0.9mm) | 5.0cm | Pink | 50 |
| 100107 | 22 (0.7mm) | 2.5cm | Blue | 50 |

Thoracic
Catheter
Accessories

Thoracic Catheter Accessories



ARGYLE EZ-FLO Stopcocks

- Blue tint polycarbonate construction.
- Low torque, smooth turning control.
- Unique threaded female luer connection.
- Sterile.

| Code 8888- | Description | No. in Carton |
|---------------|--|------------------|
| 173203 | 3-way with injection site and 51cm extension tube | 50 |
| 173500 | 3-way only | 50 |
| 173518 | 3-way with port covers | 50 |
| 174201 | 4-way with 51cm extension tube and port covers | 50 |

| Code | Product | Catalogue | Page |
|-------------|--|-----------|------|
| 1180-562009 | ARGYLE D/Lumen Thoracic Trocar Catheter, Clear PVC, 12Fr, 23cm | CT | 12 |
| 1180-562017 | ARGYLE D/Lumen Thoracic Trocar Catheter, Clear PVC, 16Fr, 23cm | CT | 12 |
| 1180-562025 | ARGYLE D/Lumen Thoracic Trocar Catheter, Clear PVC, 20Fr, 25cm | CT | 12 |
| 1180-562033 | ARGYLE D/Lumen Thoracic Trocar Catheter, Clear PVC, 24Fr, 41cm | CT | 12 |
| 1180-562041 | ARGYLE D/Lumen Thoracic Trocar Catheter, Clear PVC, 28Fr, 41cm | CT | 12 |
| 1180-571570 | THORA-SEAL II Chest Drainage Unit | CT | 1 |
| 1814-713105 | THORA-SEAL I Bottle only | CT | 1 |
| 1814-713127 | THORA-SEAL I Carrying Frame for 2 bottles | CT | 1 |
| 8884-571368 | AQUA-SEAL Evacuating Autotransfusion Refill Bag | CT | 6 |
| 8884-713100 | THORA-SEAL I Chest Drainage Unit | CT | 1 |
| 8884-713176 | THORA-Seal III Autotransfusion Chest Drainage Unit | CT | 2 |
| 8884-713184 | ARGYLE Autotransfusion receptacle for collection & reinfusion of blood | CT | 6 |
| 8884-713192 | THORA-SEAL III Autotransfusion Chest Drainage Unit w/Acc. Unit | CT | 3 |
| 8884-713308 | THORA-SEAL III Chest Drainage Unit | CT | 2 |
| 8884-713900 | Replacement Collection Chamber | CT | 3 |
| 8884-714100 | ARGYLE Tubing Extension Set, 90cm – for use w/single thoracic catheter | CT | 12 |
| 8884-714200 | ARGYLE Tubing Extension Set, 90cm – for use w/dual thoracic catheter | CT | 12 |
| 8884-722100 | CLINICAL PRODUCTS 2000ml Bottle only | CT | 10 |
| 8884-722400 | CLINICAL PRODUCTS Floor Rack | CT | 10 |
| 8884-726100 | CLINICAL PRODUCTS One Bottle Chest Drainage System, Standard | CT | 10 |
| 8884-726200 | CLINICAL PRODUCTS Two Bottle Chest Drainage System, Standard | CT | 10 |
| 8884-726400 | CLINICAL PRODUCTS Two Bottle Chest Drainage System, D/Drainage | CT | 10 |
| 8884-726500 | CLINICAL PRODUCTS One Bottle Chest Drainage System, D/Drainage | CT | 10 |
| 8884-728100 | CLINICAL PRODUCTS One Bottle Chest Drainage System, Infant | CT | 10 |
| 8884-729300 | CLINICAL PRODUCTS Fixed One Bottle Tubing Set | CT | 10 |
| 8884-729500 | CLINICAL PRODUCTS Adjustable One Bottle Tubing Set | CT | 10 |
| 8884-729600 | CLINICAL PRODUCTS Adjustable Two Bottle Tubing Set | CT | 10 |
| 8886-822726 | ARGYLE CARDIOPUNCH Disposable Aortic Punch, 2.6mm | CT | 24 |
| 8886-822740 | ARGYLE CARDIOPUNCH Disposable Aortic Punch, 4.0mm | CT | 24 |
| 8886-822745 | ARGYLE CARDIOPUNCH Disposable Aortic Punch, 4.5mm | CT | 24 |
| 8886-822748 | ARGYLE CARDIOPUNCH Disposable Aortic Punch, 4.8mm | CT | 24 |
| 8886-822750 | ARGYLE CARDIOPUNCH Disposable Aortic Punch, 5.0mm | CT | 24 |
| 8888-100107 | ARGYLE MEDICUT Intravenous Cannula, 22ga, 2.5cm | CT | 25 |
| 8888-100206 | ARGYLE MEDICUT Intravenous Cannula, 12ga, 5.0cm | CT | 25 |
| 8888-100214 | ARGYLE MEDICUT Intravenous Cannula, 14ga, 5.0cm | CT | 25 |
| 8888-100222 | ARGYLE MEDICUT Intravenous Cannula, 16ga, 5.0cm | CT | 25 |
| 8888-100230 | ARGYLE MEDICUT Intravenous Cannula, 18ga, 5.0cm | CT | 25 |
| 8888-100248 | ARGYLE MEDICUT Intravenous Cannula, 20ga, 5.0cm | CT | 25 |
| 8888-173203 | ARGYLE EZ-FLO Stopcock, 3-way w/injection site & 51cm ext. tube | CT | 26 |
| 8888-173500 | ARGYLE EZ-FLO Stopcock, 3-way only | CT | 26 |
| 8888-173518 | ARGYLE EZ-FLO Stopcock, 3-way w/covers | CT | 26 |
| 8888-173534 | ARGYLE EZ-FLO Stopcock, 3-way w/84cm ext. tube & covers | CT | 26 |
| 8888-174201 | ARGYLE EZ-FLO Stopcock, 4-way w/51cm ext. tube & covers | CT | 26 |
| 8888-513002 | ARGYLE Penrose Tubing, Latex, 6mm diam., 30cm | CT | 24 |
| 8888-513200 | ARGYLE Penrose Tubing, Latex, 8mm diam., 30cm | CT | 24 |
| 8888-513408 | ARGYLE Penrose Tubing, Latex, 10mm diam., 30cm | CT | 24 |
| 8888-513606 | ARGYLE Penrose Tubing, Latex, 13mm diam., 30cm | CT | 24 |
| 8888-513804 | ARGYLE Penrose Tubing, Latex, 16mm diam., 30cm | CT | 24 |
| 8888-514000 | ARGYLE Penrose Tubing, Latex, 19mm diam., 30cm | CT | 24 |
| 8888-514208 | ARGYLE Penrose Tubing, Latex, 22mm diam., 30cm | CT | 24 |
| 8888-514406 | ARGYLE Penrose Tubing, Latex, 25mm diam., 30cm | CT | 24 |
| 8888-514604 | ARGYLE Penrose Tubing, Latex, 6mm diam., 46cm | CT | 24 |

Index by
Code Number

| Code | Product | Catalogue | Page |
|-------------|---|-----------|------|
| 8888-514802 | ARGYLE Penrose Tubing, Latex, 8mm diam., 46cm | CT | 24 |
| 8888-515007 | ARGYLE Penrose Tubing, Latex, 10mm diam., 46cm | CT | 24 |
| 8888-515205 | ARGYLE Penrose Tubing, Latex, 13mm diam., 46cm | CT | 24 |
| 8888-515403 | ARGYLE Penrose Tubing, Latex, 16mm diam., 46cm | CT | 24 |
| 8888-515601 | ARGYLE Penrose Tubing, Latex, 19mm diam., 46cm | CT | 24 |
| 8888-515809 | ARGYLE Penrose Tubing, Latex, 25mm diam., 46cm | CT | 24 |
| 8888-530519 | ARGYLE SARATOGA Sump Drain, Clear PVC, 16Fr, 30cm | CT | 25 |
| 8888-530527 | ARGYLE SARATOGA Sump Drain, Clear PVC, 20Fr, 30cm | CT | 25 |
| 8888-530535 | ARGYLE SARATOGA Sump Drain, Clear PVC, 28Fr, 30cm | CT | 25 |
| 8888-530543 | ARGYLE SARATOGA Sump Drain, Clear PVC, 24Fr, 30cm | CT | 25 |
| 8888-530550 | ARGYLE SARATOGA Sump Drain, Clear PVC, 34Fr, 30cm | CT | 25 |
| 8888-530568 | ARGYLE SARATOGA Sump Drain, Clear PVC, 0Fr, 30cm | CT | 25 |
| 8888-531012 | ARGYLE SARATOGA Sump Drain, Clear PVC, 16Fr, 51cm | CT | 25 |
| 8888-531020 | ARGYLE SARATOGA Sump Drain, Clear PVC, 20Fr, 51cm | CT | 25 |
| 8888-531038 | ARGYLE SARATOGA Sump Drain, Clear PVC, 24Fr, 51cm | CT | 25 |
| 8888-531046 | ARGYLE SARATOGA Sump Drain, Clear PVC, 28Fr, 51cm | CT | 25 |
| 8888-531053 | ARGYLE SARATOGA Sump Drain, Clear PVC, 34Fr, 51cm | CT | 25 |
| 8888-531061 | ARGYLE SARATOGA Sump Drain, Clear PVC, 40Fr, 51cm | CT | 25 |
| 8888-551010 | ARGYLE Gott Aortic Aneurysm Shunt, Clear PVC, 66cm | CT | 22 |
| 8888-551036 | ARGYLE Gott Aortic Aneurysm Shunt, Clear PVC, 90cm | CT | 22 |
| 8888-552018 | ARGYLE Hepatic Venous Shunt, Clear PVC, Adult | CT | 22 |
| 8888-560805 | ARGYLE Trocar Catheter, Clear PVC, 8Fr, 23cm | CT | 12 |
| 8888-561019 | ARGYLE Trocar Catheter, Clear PVC, 10Fr, 23cm | CT | 12 |
| 8888-561027 | ARGYLE Trocar Catheter, Clear PVC, 12Fr, 23cm | CT | 12 |
| 8888-561035 | ARGYLE Trocar Catheter, Clear PVC, 16Fr, 25cm | CT | 12 |
| 8888-561043 | ARGYLE Trocar Catheter, Clear PVC, 20Fr, 41cm | CT | 12 |
| 8888-561050 | ARGYLE Trocar Catheter, Clear PVC, 24Fr, 41cm | CT | 12 |
| 8888-561068 | ARGYLE Trocar Catheter, Clear PVC, 28Fr, 41cm | CT | 12 |
| 8888-561076 | ARGYLE Trocar Catheter, Clear PVC, 32Fr, 41cm | CT | 12 |
| 8888-561514 | ARGYLE Trocar Replacement Catheter, Clear PVC, 10Fr, 23cm | CT | 12 |
| 8888-561522 | ARGYLE Trocar Replacement Catheter, Clear PVC, 12Fr, 23cm | CT | 12 |
| 8888-561530 | ARGYLE Trocar Replacement Catheter, Clear PVC, 16Fr, 23cm | CT | 12 |
| 8888-561548 | ARGYLE Trocar Replacement Catheter, Clear PVC, 20Fr, 25cm | CT | 12 |
| 8888-561555 | ARGYLE Trocar Replacement Catheter, Clear PVC, 24Fr, 41cm | CT | 12 |
| 8888-561563 | ARGYLE Trocar Replacement Catheter, Clear PVC, 28Fr, 41cm | CT | 12 |
| 8888-561571 | ARGYLE Trocar Replacement Catheter, Clear PVC, 32Fr, 41cm | CT | 12 |
| 8888-565028 | ARGYLE Trocar Catheter Kit, 12Fr, 23cm | CT | 15 |
| 8888-565036 | ARGYLE Trocar Catheter Kit, 20Fr, 25cm | CT | 15 |
| 8888-566067 | ARGYLE TURKEL Thoracentesis Set, 8Fr | CT | 16 |
| 8888-570507 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 12Fr, 51cm | CT | 11 |
| 8888-570515 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 16Fr, 51cm | CT | 11 |
| 8888-570523 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 20Fr, 51cm | CT | 11 |
| 8888-570531 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 24Fr, 51cm | CT | 11 |
| 8888-570549 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 28Fr, 51cm | CT | 11 |
| 8888-570556 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 32Fr, 51cm | CT | 11 |
| 8888-570564 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 36Fr, 51cm | CT | 11 |
| 8888-570572 | ARGYLE Thoracic Catheter, Clear PVC, Straight, 40Fr, 51cm | CT | 11 |
| 8888-571018 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 16Fr, 51cm | CT | 11 |
| 8888-571026 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 20Fr, 51cm | CT | 11 |
| 8888-571034 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 24Fr, 51cm | CT | 11 |
| 8888-571042 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 28Fr, 51cm | CT | 11 |
| 8888-571059 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 32Fr, 51cm | CT | 11 |
| 8888-571067 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 36Fr, 51cm | CT | 11 |
| 8888-571075 | ARGYLE Thoracic Catheter, Clear PVC, Right Angle, 40Fr, 51cm | CT | 11 |
| 8888-571299 | AQUA-SEAL Chest Drainage Unit | CT | 4 |
| 8888-571315 | AQUA-SEAL Autotransfusion Chest Drainage Unit | CT | 5 |
| 8888-571323 | AQUA-SEAL Autotransfusion Chest Drainage Unit w/Acc. Unit | CT | 5 |
| 8888-571364 | AQUA-SEAL All Purpose Autotransfusion Chest Drainage Unit | CT | 6 |
| 8888-571406 | AQUA-SEAL Dual Drain Chest Drainage Unit | CT | 4 |
| 8888-571463 | SENTINEL SEAL Autotransfusion Chest Drainage Unit w/Acc. Unit | CT | 8 |
| 8888-571489 | SENTINEL SEAL Autotransfusion Chest Drainage Unit | CT | 8 |

Index by Code Number

| Code | Product | Catalogue | Page |
|-------------|--|-----------|------|
| 8888-571513 | SENTINEL SEAL Dual Drain Chest Drainage Unit | CT | 7 |
| 8888-571547 | DOUBLE SEAL Chest Drainage Unit | CT | 9 |
| 8888-571562 | SENTINEL SEAL Chest Drainage Unit | CT | 7 |
| 8888-571588 | SINGLE SEAL Chest Drainage Unit, 2000ml, Two Bottle Unit | CT | 1 |
| | | | |
| 8888-572537 | ARGYLE Thoracic Catheter, Silicone, Straight, 20Fr, 51cm | CT | 11 |
| 8888-572545 | ARGYLE Thoracic Catheter, Silicone, Straight, 24Fr, 51cm | CT | 11 |
| 8888-572552 | ARGYLE Thoracic Catheter, Silicone, Straight, 28Fr, 51cm | CT | 11 |
| 8888-572560 | ARGYLE Thoracic Catheter, Silicone, Straight, 32Fr, 51cm | CT | 11 |
| 8888-572578 | ARGYLE Thoracic Catheter, Silicone, Straight, 36Fr, 51cm | CT | 11 |
| 8888-572586 | ARGYLE Thoracic Catheter, Silicone, Straight, 40Fr, 51cm | CT | 11 |
| | | | |
| 8888-573030 | ARGYLE Thoracic Catheter, Silicone, Right Angle, 24Fr, 51cm | CT | 11 |
| 8888-573048 | ARGYLE Thoracic Catheter, Silicone, Right Angle, 28Fr, 51cm | CT | 11 |
| 8888-573055 | ARGYLE Thoracic Catheter, Silicone, Right Angle, 32Fr, 51cm | CT | 11 |
| 8888-573063 | ARGYLE Thoracic Catheter, Silicone, Right Angle, 36Fr, 51cm | CT | 11 |
| 8888-573550 | ARGYLE Mediastinal Cath., Silicone, 28Fr, 51cm | CT | 12 |
| 8888-573568 | ARGYLE Mediastinal Cath., Silicone, 34Fr, 51cm | CT | 12 |
| | | | |
| 8888-577775 | ARGYLE Carotid Artery Shunt Kit, 4 Tubes: 8/10/12/14Fr, 15cm | CT | 22 |
| 8888-584912 | ARGYLE Parallel "Y" Connector, 9.5 x 9.5 x 13mm | CT | 20 |
| 8888-585000 | ARGYLE Vascular Tourniquet Kit, 13cm | CT | 23 |
| 8888-585208 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Regular, 38.0cm | CT | 23 |
| 8888-585216 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Regular, 38.0cm | CT | 23 |
| 8888-585224 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Regular, 38.0cm | CT | 23 |
| 8888-585232 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Regular, 38.0cm | CT | 23 |
| 8888-585240 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Small, 35.5cm | CT | 23 |
| 8888-585257 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Small, 35.5cm | CT | 23 |
| 8888-585265 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Small, 35.5cm | CT | 23 |
| 8888-585273 | ARGYLE SENTINEL LOOPS Vascular Tape, Silicone, Small, 35.5cm | CT | 23 |
| 8888-590018 | ARGYLE Ferguson Left Vent Catheter, 18Fr, 51cm | CT | 19 |
| 8888-590026 | ARGYLE Ferguson Left Vent Catheter, 24Fr, 51cm | CT | 19 |
| 8888-590208 | ARGYLE Left Ventricular Sump Vent Catheter, 14Fr, 51cm | CT | 19 |
| 8888-590216 | ARGYLE Left Ventricular Sump Vent Catheter, 18Fr, 51cm | CT | 19 |
| 8888-590224 | ARGYLE Left Ventricular Sump Vent Catheter, 20Fr, 51cm | CT | 19 |
| 8888-591008 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 10Fr | CT | 20 |
| 8888-591016 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 10Fr | CT | 20 |
| 8888-591024 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 14Fr | CT | 20 |
| 8888-591032 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 14Fr | CT | 20 |
| 8888-591040 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 18Fr | CT | 20 |
| 8888-591057 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 18Fr | CT | 20 |
| 8888-591065 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 21Fr | CT | 20 |
| 8888-591073 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 21Fr | CT | 20 |
| 8888-591081 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 24Fr | CT | 20 |
| 8888-591099 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 24Fr | CT | 20 |
| 8888-591123 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 12Fr | CT | 20 |
| 8888-591131 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 12Fr | CT | 20 |
| 8888-591149 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Angled, 16Fr | CT | 20 |
| 8888-591156 | ARGYLE THI Aortic Perfusion Cannula, Clear PVC, Straight, 16Fr | CT | 20 |
| 8888-593012 | ARGYLE THI Perfusion Adapter, Adult | CT | 20 |
| 8888-593038 | ARGYLE THI Perfusion Adapter, Paediatric | CT | 20 |
| 9999-571548 | MINIDRAIN I Pump, 6L | CT | 21 |
| 9999-571558 | MINIDRAIN II Pump, 11L | CT | 21 |

Index by
Code Number

Altitude*

Thoracic drainage system with three chambers and dry suction adjustment

* Trademark © 2002 Tyco Healthcare Group LP. All rights reserved.



Altitude* - Product specifications

| Ref | Description | Collection capacity | Qty per package |
|-------------|--|---------------------|-----------------|
| 8888-571370 | Altitude* pre-set for auto-transfusion | 2500 ml | 5 |

* Trademark © 2004 Tyco Healthcare Group LP. All rights reserved

tyco
Healthcare

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT Hampshire
PO 13 0AS
UK
Tel: +44 (0) 1329 224226
Fax: +44 (0) 1329 224334

Tyco Healthcare Regional
Export
154 Fareham Road
GOSPORT Hampshire
PO 13 0AS
UK
Tel: +44 (0) 1329 224258
Fax: +44 (0) 1329 224086

Tyco Healthcare Nederland BV
Hogeweg 105
NL-5301 LL ZALTBOMMEL
The NETHERLANDS
Tel: +31 418 57 66 00
Fax: +31 418 57 67 91
TycoHealthNed@TycoHealth.com
www.tycohealthcare.nl

Tyco Healthcare
Birchwood Court Montrose Road
P.O. Box 8108, Halfway House
1685 MIDRAND
SOUTH AFRICA
Tel: +27 (0) 11 265 2000
Fax: +27 (0) 11 265 2050

Tyco Healthcare Belgium N.V.-S.A.
Generaal de Wittelaan 9/5
B-2800 MECHELEN
BELGIUM
Tel: +32 (0) 15 29 44 50
Fax: +32 (0) 15 29 44 55
tycohealthcare.belgium@tyco-healthcare.com

Tyco Healthcare Norden AB
Gårdsvägen 18, Box 710
SE-169 27 Solna
SWEDEN
Tel: +46-8 585 605 00
Kundservice +46-8 585 605 05
Fax: +46-8 585 605 06

Tyco Healthcare Danmark
Langebrogade 6E, 4. th.
DK-1411 København K
DENMARK
Tel: +45-70 27 52 50
Kundeservice +45-70 27 52 50
Fax: +45-70 27 51 50

Tyco Healthcare Norge
Løkketangen 20
N-1337 Sandvika
NORWAY
Tel: +47-67 80 55 10
Kundeservice +47-67 80 55 18
Fax: +47-67 80 55 11

Tyco Healthcare Finland
Köydenpunojankatu 2D, 6 krs.
PL 340
FIN-00181 Helsinki
FINLAND
Puh +358-(0)9-6226 840
Asiakaspalvelu
+358-(0)9-6226 8418
Fax: +358-(0)9-6226 8460