



TELFA[™] AMD[™]

Non - adherent dressing



- Preventative protection against nosocomial wound infections.
- AMD provides broad spectrum effectiveness whilst being gentle to healthy cells.
- Helps to provide a Low Cost Solution to Infection Control.

TELFA AMD is the first product to be launched within the AMD range - more will follow!



TELFA™ AMD™ Dressings



AMD helps to prevent colonisation - Prevention is better than Treatment

What is AMD?

The antimicrobial component PHMB (Polyhexamethylene Biguanide) resists bacterial colonization with Gram-positive, Gram-negative and Fungi/Yeast microorganisms within and through the dressing, helping to reduce bacterial colonisation.

Reduce infection - Reduce Costs

Infections increase your costs significantly - AMD helps to provide low cost prevention against nosocomial infection, helping to reduce the length of hospital stay.

Unique & effective

- Broad spectrum effectiveness
- No known resistance
- Works within and through the dressing
- Gentle to healthy cells

Main Indications

- Surgical Incisions
- Suture sites
- Abrasion
- Lacerations
- Dry to light exudating wounds
- I.V. lines
- Arterial stasis ulcers
- Helps to prevent wound infections

TELFA AMD Island Dressings have demonstrable in-vitro effectiveness against:

Staphylococcus aureus
Enterococcus faecalis
Acinetobacter baumannii
Candida albicans
Escherichia coli
Pseudomonas aeruginosa
Staphylococcus epidermidis
Enterobacter cloacae
Klebsiella pneumoniae
Proteus mirabilis
Serratia marcescens
Staphylococcus coagulase

Test data on file at Kendall

Ordering Information

TELFA AMD Pads Sterile

| Code | Description | Per Tray | Ship Case |
|------|-------------|----------|-----------|
| 7662 | 7.5 x 10 cm | 50 | 900 |
| 7663 | 7.5 x 20 cm | 50 | 600 |

TELFA AMD Island Dressings

| Code | Description | Per Tray | Ship Case |
|------|--------------|----------|-----------|
| 7665 | 10 x 12.5 cm | 25 | 200 |
| 7666 | 10 x 20 cm | 25 | 100 |
| 7667 | 10 x 25.5 cm | 25 | 100 |
| 7668 | 10 x 35 cm | 25 | 50 |

Any further Questions? contact me directly at:
Amadeus@emea.tycohealthcare.com



tyco
Healthcare

Kendall

TELFA™ and AMD™ are trademarks of Tyco Healthcare Group LP. © 2004 Tyco Healthcare Group LP. All rights reserved.

Tyco Healthcare UK Limited
154 Fareham Road
GOSPORT - HAMPSHIRE
UK - PO13 0AS
Tel: +44 (0) 13 29 22 42 26
Fax: +44 (0) 13 29 22 43 34